

2-Propanol  
HPLC  
For use in Liquid Chromatography (HPLC & UHPLC) &  
Spectrophotometry

**LabGenius**

Material No: 1.001001  
Batch No.: G0163311  
Manufactured Date: 2024/07/01  
Retest Date: 2029/07/01

## Certificate of Analysis

| Test   | Specification | Result |
|--|---------------|--------|
| Ultraviolet Absorbance (1.00-cm cell vs. water) – 350 nm         | <= 0.01       | < 0.01 |
| Ultraviolet Absorbance (1.00-cm cell vs. water) – 280 nm         | <= 0.01       | < 0.01 |
| Ultraviolet Absorbance (1.00-cm cell vs. water) – 254 nm         | <= 0.020      | 0.010  |
| Ultraviolet Absorbance (1.00-cm cell vs. water) – 225 nm         | <= 0.16       | 0.10   |
| Ultraviolet Absorbance (1.00-cm cell vs. water) – UV Cut-off, nm | <= 205        | 204    |
| Assay (CH <sub>3</sub> CHOHCH <sub>3</sub> )                     | >= 99.7 %     | 99.9   |
| Residue after Evaporation  | <= 2.0 ppm    | 0.9    |
| Water (by KF, coulometric)                                       | <= 0.03 %     | 0.01   |

For Laboratory, Research or Manufacturing Use

Country of Origin: CN  
Packaging Site: Paris Mfg Ctr & DC